



MUC16/CA125 Antibody

Product Code	CSB-PA899266
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WXI7
Immunogen	Synthesized peptide derived from internal of Human MUC16.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,IF;IF:1:100-1:500
Form	Liquid
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MUC16

Image



Immunofluorescence analysis of HepG2 cells, using MUC16 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.