




# MUC16/CA125 Antibody

<b>Product Code</b>	CSB-PA899266
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q8WXI7
<b>Immunogen</b>	Synthesized peptide derived from internal of Human MUC16.
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA,IF;IF:1:100-1:500
<b>Form</b>	Liquid
<b>Purification Method</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Clonality</b>	Polyclonal
<b>Product Type</b>	Polyclonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	MUC16

## Image



Immunofluorescence analysis of HepG2 cells, using MUC16 antibody.