




TFAP2C Antibody

Product Code	CSB-PA036945
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92754
Immunogen	Synthesized peptide derived from internal of Human AP2C.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total AP2C protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	AP2-gamma; Activating enhancer-binding protein 2 gamma; Transcription factor ERF-1;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TFAP2C

Image



Western blot analysis of extracts from HepG2 cells and HUVEC cells, using AP2C antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.