




CELF1 Antibody

Product Code	CSB-PA991834
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92879
Immunogen	Synthesized peptide derived from internal of Human CELF-1.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total CELF-1 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	CELF-1; Bruno-like protein 2; RNA-binding protein BRUNOL-2; CUG triplet repeat RNA-binding protein 1; CUG-BP1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CELF1

Image



Western blot analysis of extracts from HeLa cells and HepG2 cells, using CELF-1 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.