




PHKB Antibody

Product Code	CSB-PA901833
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q93100
Immunogen	Synthesized peptide derived from internal of Human KPBB.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total KPBB protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Phosphorylase b kinase regulatory subunit beta; Phosphorylase kinase subunit beta; PHKB;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PHKB

Image



Western blot analysis of extracts from K562 cells, using KPBB antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.