




C18orf8 Antibody

Product Code	CSB-PA113796
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96DM3
Immunogen	Synthesized peptide derived from C-terminal of Human MIC1.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total MIC1 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Uncharacterized protein C18orf8; Colon cancer-associated protein Mic1; Mic-1; C18orf8;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	RMC1

Image



Western blot analysis of extracts from Jurkat cells and HepG2 cells, using MIC1 antibody.