




HES6 Antibody

Product Code	CSB-PA788871
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96HZ4
Immunogen	Synthesized peptide derived from N-terminal of Human HES6.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total HES6.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Transcription cofactor HES-6; Hairy and enhancer of split 6; C-HAIRY1; HES6;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	HES6

Image



Western blot analysis of extracts from HeLa cells, using HES6 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.