




FOXB1/FOXB2 Antibody

Product Code	CSB-PA027895
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99853
Immunogen	Synthesized peptide derived from internal of Human FOXB1/2.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous levels of total FOXB1/2 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Forkhead box protein B1/2; Transcription factor FKH-5; FOXB1/2; FKH5;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	FOXB1/FOXB2

Image



Western blot analysis of extracts from 293 cells, using FOXB1/2 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.