




ADPGK Antibody

Product Code	CSB-PA272036
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BRR6
Immunogen	Synthesized peptide derived from internal of Human ADPGK.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total ADPGK protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	ADP-dependent glucokinase; EC 2.7.1.147; ADPGK; ADP-GK; RbBP-35
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ADPGK

Image



Western blot analysis of extracts from Jurkat cells, using ADPGK antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.