




MKNK1 Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA167602 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q9BUB5 |
| Immunogen | Synthesized peptide derived from internal of Human MKNK1. |
| Raised In | Rabbit |
| Species Reactivity | Human,Mouse |
| Specificity | The antibody detects endogenous levels of total MKNK1 protein. |
| Tested Applications | ELISA,WB;WB:1:500-1:3000 |
| Form | Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Purification Method | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Clonality | Polyclonal |
| Alias | MAP kinase-interacting serine/threonine-protein kinase 1; EC 2.7.11.1; MAP kinase signal-integrating kinase 1; Mnk1; MKNK1 |
| Product Type | Polyclonal Antibody |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | MKNK1 |

Image



Western blot analysis of extracts from HuvEc cells, using MKNK1 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.