




NAA15 Antibody

Product Code	CSB-PA058433
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BXJ9
Immunogen	Synthesized peptide derived from internal of Human NARG1.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total NARG1 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	N-terminal acetyltransferase; Protein tubedown-1; Tbdn100; Gastric cancer antigen Ga19;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	NAA15

Image



Western blot analysis of extracts from HUVEC cells, using NARG1 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.