




IFIH1 Antibody

Product Code	CSB-PA255539
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BYX4
Immunogen	Synthesized peptide derived from C-terminal of Human IFIH1.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total IFIH1 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Interferon-induced helicase C domain-containing protein 1; EC 3.6.1.-; Interferon-induced with helicase C domain protein 1; Helicase with 2 CARD domains; Helicard
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	IFIH1

Image



Western blot analysis of extracts from Jurkat cells, using IFIH1 antibody.