




PKIB Antibody

Product Code	CSB-PA179075
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9C010
Immunogen	Synthesized peptide derived from internal of Human IPKB.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous levels of total IPKB protein.
Tested Applications	ELISA, IHC, IF; IHC: 1:50-1:100, IF: 1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	PKI-beta; cAMP-dependent protein kinase inhibitor beta; PKIB; PRKACN2;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PKIB

Image



Immunohistochemistry analysis of paraffin-embedded human brain tissue, using IPKB antibody.



Immunofluorescence analysis of MCF-7 cells, using IPKB antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.