





C2orf40 Antibody

Product Code	CSB-PA067339
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H1Z8
Immunogen	Synthesized peptide derived from internal of Human ECRG4.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous levels of total ECRG4 protein.
Tested Applications	ELISA, IHC, IF; IHC: 1:50-1:100, IF: 1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Esophageal cancer-related gene 4 protein [Precursor]; C2orf40; ECRG4; UNQ761/PRO1508;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	C2orf40

Image



Immunohistochemistry analysis of paraffin-embedded human testis tissue, using ECRG4 antibody.



Immunofluorescence analysis of A549 cells, using ECRG4 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.