





MUC13 Antibody

Product Code	CSB-PA276815
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H3R2
Immunogen	Synthesized peptide derived from internal of Human MUC13.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total MUC13 protein.
Tested Applications	ELISA,WB,IHC,IF;WB:1:500-1:3000,IHC:1:50-1:100,IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Mucin-13 [Precursor]; Down-regulated in colon cancer 1; MUC13; DRCC1; RECC
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MUC13

Image



Western blot analysis of extracts from COLO cells, 293 cells and HeLa cells, using MUC13 antibody.



Immunohistochemistry analysis of paraffin-embedded human colon carcinoma tissue using MUC13 antibody.



Immunofluorescence analysis of HepG2 cells, using MUC13 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.