




FKRP Antibody

Product Code	CSB-PA035204
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H9S5
Immunogen	Synthesized peptide derived from N-terminal of Human FKRP.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total FKRP protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Fukutin-related protein; EC 2.-.-.-; FKRP;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	FKRP

Image



Western blot analysis of extracts from Jurkat cells and COLO205 cells, using FKRP antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.