




FKRP Antibody

Product Code	CSB-PA035204
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H9S5
Immunogen	Synthesized peptide derived from N-terminal of Human FKRP.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total FKRP protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Fukutin-related protein; EC 2.-.-.-; FKRP;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	FKRP

Image



Western blot analysis of extracts from Jurkat cells and COLO205 cells, using FKRP antibody.