





ZNF287 Antibody

Product Code	CSB-PA555215
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HBT7
Immunogen	Synthesized peptide derived from internal of Human ZNF287.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total ZNF287 protein.
Tested Applications	ELISA, WB, IHC; WB: 1:500-1:3000, IHC: 1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Zinc finger protein with KRAB and SCAN domains 13; Zinc finger protein 287; ZNF287; ZKSCAN13;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ZNF287

Image



Western blot analysis of extracts from NIH-3T3 cells, using ZNF287 antibody.



Immunohistochemistry analysis of paraffin-embedded human brain tissue using ZNF287 antibody.