




ARNT2 Antibody

Product Code	CSB-PA133307
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HBZ2
Immunogen	Synthesized peptide derived from internal of Human ARNT2.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total ARNT2 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	ARNT protein 2; Aryl hydrocarbon receptor nuclear translocator 2; ARNT2; KIAA0307;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ARNT2

Image



Western blot analysis of extracts from Jurkat cells and 293 cells, using ARNT2 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.