





AKIP1 Antibody

Product Code	CSB-PA956447
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NQ31
Immunogen	Synthesized peptide derived from internal of Human BCA3.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total BCA3 protein.
Tested Applications	ELISA, WB, IHC; WB: 1:500-1:3000, IHC: 1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Breast cancer-associated gene 3 protein; PKA-interacting protein; AKIP1;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	AKIP1

Image



Western blot analysis of extracts from Jurkat cells, using BCA3 antibody.



Immunohistochemistry analysis of paraffin-embedded human heart tissue using BCA3 antibody.