




ETNK2 Antibody

Product Code	CSB-PA029059
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NVF9
Immunogen	Synthesized peptide derived from internal of Human EK12.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total EK12 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Ethanolamine kinase 2; EC 2.7.1.82; EK1 2; Ethanolamine kinase-like protein; ETNK2
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ETNK2

Image



Western blot analysis of extracts from 293 cells, using EK12 antibody.