




DNAJC17 Antibody

Product Code	CSB-PA053434
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NVM6
Immunogen	Synthesized peptide derived from N-terminal of Human DNAJC17.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total DNAJC17 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	DnaJ homolog subfamily C member 17;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	DNAJC17

Image



Western blot analysis of extracts from Jurkat cells, HepG2 cells and K562 cells, using DNAJC17 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.