




AATF Antibody

Product Code	CSB-PA970851
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NY61
Immunogen	Synthesized peptide derived from N-terminal of Human AATF.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total AATF protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	CHE-1; CHE1; DED; apoptosis antagonizing transcription factor;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	AATF

Image



Western blot analysis of extracts from HeLa cells, HepG2 cells, COLO205 cells and HUVEC cells, using AATF antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.