




ERAP1 Antibody

Product Code	CSB-PA987708
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NZ08
Immunogen	Synthesized peptide derived from internal of Human ARTS-1.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total ARTS-1 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	A- LAP; APPILS; ARTS-1; Adipocyte-derived leucine aminopeptidase precursor; Aminopeptidase PILS
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ERAP1

Image



Western blot analysis of extracts from HepG2 cells, Jurkat cells and HUVEC cells, using ARTS-1 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.