




EIF2AK3 Antibody

Product Code	CSB-PA590642
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NZJ5
Immunogen	Synthesized peptide derived from internal of Human E2AK3.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total E2AK3 protein.
Tested Applications	ELISA,WB,IF;WB:1:500-1:3000,IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	E2AK3; EC 2.7.11.1; EIF2AK3; Eukaryotic translation initiation factor 2-alpha kinase 3 precursor; HsPEK
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	EIF2AK3

Image



Western blot analysis of extracts from HeLa cells and HepG2 cells, using E2AK3 antibody.



Immunofluorescence analysis of HepG2 cells, using E2AK3 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.