




HMG20B Antibody

Product Code	CSB-PA788036
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9P0W2
Immunogen	Synthesized peptide derived from N-terminal of Human HMG20B.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total HMG20B protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	HMGX1; high-mobility group 20A;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	HMG20B

Image



Western blot analysis of extracts from Jurkat cells, using HMG20B antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.