




DBF4 Antibody

Product Code	CSB-PA113425
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UBU7
Immunogen	Synthesized peptide derived from N-terminal of Human DBF4.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous levels of total DBF4 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Protein DBF4 homolog A; Activator of S phase kinase; DBF4-type zinc finger-containing protein 1; Chiffon homolog A; DBF4
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	DBF4

Image



Western blot analysis of extracts from NIH-3T3 cells treated with H₂O₂ (100uM, 30mins) and COS-7 cells treated with PMA (125mg/ml, 30mins), using DBF4 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.