




NXF1 Antibody

Product Code	CSB-PA260303
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UBU9
Immunogen	Synthesized peptide derived from N-terminal of Human NXF1.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total NXF1 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	MEX67; TAP; Tip-associated protein; mRNA export factor TAP; nuclear RNA export factor 1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	NXF1

Image



Western blot analysis of extracts from HepG2 cells, Jurkat cells and HeLa cells, using NXF1 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.