




## **POLL** Antibody

<b>Product Code</b>	CSB-PA100387
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UGP5
Immunogen	Synthesized peptide derived from internal of Human DNA Polymerase $\boldsymbol{\lambda}$
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total DNA Polymerase $\boldsymbol{\lambda}$ protein.
<b>Tested Applications</b>	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	DNA polymerase lambda; EC 2.7.7.7; EC 4.2.99; Pol Lambda; DNA polymerase kappa
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	POLL





Western blot analysis of extracts from 293 cells, using DNA Polymerase 1>> antibody.