




POLL Antibody

Product Code	CSB-PA100387
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UGP5
Immunogen	Synthesized peptide derived from internal of Human DNA Polymerase λ
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total DNA Polymerase λ protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	DNA polymerase lambda; EC 2.7.7.7; EC 4.2.99.-; Pol Lambda; DNA polymerase kappa
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	POLL

Image



Western blot analysis of extracts from 293 cells, using DNA Polymerase λ antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.