




GTF2IRD1 Antibody

Product Code	CSB-PA911147
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UHL9
Immunogen	Synthesized peptide derived from N-terminal of Human GTF2IRD1.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total GTF2IRD1 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	CREAM1; G2D1; GT2D1; GTF2I repeat domain containing protein 1; GTF2IRD1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	GTF2IRD1

Image



Western blot analysis of extracts from 293 cells, using GTF2IRD1 antibody.