




GTF2IRD1 Antibody

Product Code	CSB-PA911147
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UHL9
Immunogen	Synthesized peptide derived from N-terminal of Human GTF2IRD1.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Specificity	The antibody detects endogenous levels of total GTF2IRD1 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	CREAM1; G2D1; GT2D1; GTF2I repeat domain containing protein 1; GTF2IRD1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	GTF2IRD1

Image



Western blot analysis of extracts from 293 cells, using GTF2IRD1 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.