




CNOT7 Antibody

Product Code	CSB-PA446150
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UIV1
Immunogen	Synthesized peptide derived from N-terminal of Human CNOT7.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous levels of total CNOT7 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	BTG1 binding factor 1; BTG1-binding factor 1; CCR4-NOT transcription complex; 7; CCR4-associated factor 1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CNOT7

Image



Western blot analysis of extracts from HUVEC cells and HeLa cells, using CNOT7 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.