




STAG3 Antibody

Product Code	CSB-PA167892
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UJ98
Immunogen	Synthesized peptide derived from C-terminal of Human STAG3.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total STAG3 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Cohesin subunit SA-3; Stromal antigen 3; Stromalin 3; SCC3 homolog 3; STAG3
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	STAG3

Image



Western blot analysis of extracts from Jurkat cells, using STAG3 antibody.