





FOXD3 Antibody

Product Code	CSB-PA298737
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UJU5
Immunogen	Synthesized peptide derived from internal of Human FOXD3.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total FOXD3 protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Forkhead box protein D3; HNF3/FH transcription factor genesis; FOXD3;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	FOXD3

Image



Western blot analysis of extracts from HUVEC cells and 293 cells, using FOXD3 antibody.



Immunohistochemistry analysis of paraffin-embedded human heart tissue using FOXD3 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.