




ELF5 Antibody

Product Code	CSB-PA799543
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UKW6
Immunogen	Synthesized peptide derived from internal of Human ELF5.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous levels of total ELF5 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	ETS-related transcription factor Elf-5; E74-like factor 5; Epithelium-specific Ets transcription factor 2; ESE-2; Epithelium-restricted ESE-1-related Ets factor
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ELF5

Image



Western blot analysis of extracts from 293 cells, Jurkat cells and HUVEC cells, using ELF5 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.