





TPX2 Antibody

Product Code	CSB-PA560765
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9ULW0
Immunogen	Synthesized peptide derived from internal of Human DIL-2.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total DIL-2 protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	C20orf1; C20orf2; DIL-2; DIL2; Differentially expressed in lung cells 2
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TPX2

Image



Western blot analysis of extracts from Jurkat cells, using DIL-2 antibody.



Immunohistochemistry analysis of paraffin-embedded human brain tissue using DIL-2 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.