





MOK Antibody

Product Code	CSB-PA105552
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UQ07
Immunogen	Synthesized peptide derived from internal of Human MOK.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total MOK protein.
Tested Applications	ELISA, WB, IHC; WB: 1:500-1:3000, IHC: 1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	EC 2.7.11.22; MOK protein kinase; RAGE; RAGE-1; RAGE1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MOK

Image



Western blot analysis of extracts from HepG2 cells, COLO205 cells and 293 cells, using MOK antibody.



Immunohistochemistry analysis of paraffin-embedded human brain tissue using MOK antibody.