




PIKFYVE Antibody

Product Code	CSB-PA072525
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y2I7
Immunogen	Synthesized peptide derived from internal of Human PIP5K.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total PIP5K protein.
Tested Applications	ELISA,WB,IF;WB:1:500-1:3000,IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	1-phosphatidylinositol-4-phosphate 5-kinase; EC 2.7.1.68; FYV1; FYVE finger-containing phosphoinositide kinase; KIAA0981
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PIKFYVE

Image



Western blot analysis of extracts from HepG2 cells, using PIP5K antibody.



Immunofluorescence analysis of COS7 cells, using PIP5K antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.