



ZKSCAN5 Antibody

Product Code	CSB-PA221980
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y2L8
Immunogen	Synthesized peptide derived from internal of Human ZNF95.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total ZNF95 protein.
Tested Applications	ELISA,IF;IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Zinc finger protein with KRAB and SCAN domains 5; Zinc finger protein 95 homolog; Zfp-95; ZKSCAN5; KIAA1015
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ZKSCAN5

Image



Immunofluorescence analysis of A549 cells, using ZNF95 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.