




DYRK1B Antibody

Product Code	CSB-PA908230
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y463
Immunogen	Synthesized peptide derived from internal of Human DYR1B.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous levels of total DYR1B protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	DYR1B; Dual-specificity tyrosine-phosphorylation regulated kinase 1B; EC 2.7.12.1; MIRK; Mirk protein kinase
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	DYRK1B

Image



Western blot analysis of extracts from Jurkat cells, using DYR1B antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.