





ATG4B Antibody

Product Code	CSB-PA138409
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y4P1
Immunogen	Synthesized peptide derived from internal of Human ATG4B.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous levels of total ATG4B protein.
Tested Applications	ELISA, WB, IHC; WB: 1:500-1:3000, IHC: 1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Cysteine protease ATG4B; EC 3.4.22.-; Autophagy-related protein 4 homolog B; hAPG4B; Autophagin-1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ATG4B

Image



Western blot analysis of extracts from Jurkat cells, using ATG4B antibody.



Immunohistochemistry analysis of paraffin-embedded human skeletal muscle tissue using ATG4B antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.