





CDC42BPB Antibody

Product Code	CSB-PA936608
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y5S2
Immunogen	Synthesized peptide derived from internal of Human MRCKB.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total MRCKB protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	CDC42-binding protein kinase beta; CDC42BPB; EC 2.7.11.1; MRCK-beta;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CDC42BPB

Image



Western blot analysis of extracts from COLO cells, using MRCKB antibody.



Immunohistochemistry analysis of paraffin-embedded human brain tissue using MRCKB antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.