




TAF6L Antibody

Product Code	CSB-PA438797
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y6J9
Immunogen	Synthesized peptide derived from internal of Human TAF6L.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total TAF6L protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	TAF6-like RNA polymerase II p300/CBP-associated factor-associated factor 65 kDa subunit 6L; PCAF-associated factor 65 alpha; PAF65-alpha; TAF6L; PAF65A
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TAF6L

Image



Western blot analysis of extracts from 293 cells, using TAF6L antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.