




IKBKG Antibody

Product Code	CSB-PA055629
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y6K9
Immunogen	Synthesized peptide derived from internal of Human IKK- γ .
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total IKK- γ protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	FIP-3; FIP3; I-kappa-B kinase gamma; IKBKG; IKK-gamma
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	IKBKG

Image



Western blot analysis of extracts from COS-7 cells, using IKK- γ antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.