




ETV4 Antibody

Product Code	CSB-PA096475
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P43268
Immunogen	Synthesized peptide derived from internal of Human ETV4.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous levels of total ETV4 protein.
Tested Applications	ELISA, IHC; IHC: 1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	ETS translocation variant 4; Adenovirus E1A enhancer-binding protein; E1A-F; ETV4; E1AF
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ETV4

Image



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue, using ETV4 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.