




ADCY4 Antibody

Product Code	CSB-PA160118
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NFM4
Immunogen	Synthesized peptide derived from internal of Human ADCY4.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total ADCY4 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Adenylate cyclase; type IV; Adenylyl cyclase 4; ATP pyrophosphate-lyase 4; EC 4.6.1.1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ADCY4

Image



Western blot analysis of extracts from Jurkat cells, using ADCY4 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.