





ADCY6/ADCY5 Antibody

Product Code	CSB-PA108140
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43306
Immunogen	Synthesized peptide derived from internal of Human ADCY5/6.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous levels of total ADCY5/6 protein.
Tested Applications	ELISA, IHC, IF; IHC: 1:50-1:100, IF: 1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	ATP pyrophosphate-lyase; CYA6; Ca(2 ⁺)-inhibitible adenylyl cyclase; EC 4.6.1.1; KIAA0422
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ADCY6/ADCY5

Image



Immunohistochemistry analysis of paraffin-embedded human brain tissue using ADCY5/6 antibody.



Immunofluorescence analysis of HeLa cells, using ADCY5/6 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.