




ADD3 Antibody

Product Code	CSB-PA964364
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UEY8
Immunogen	Synthesized peptide derived from internal of Human ADD3.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total ADD3 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	ADDG; ADDL; adducin; gamma; adducin-like protein 70
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ADD3

Image



Western blot analysis of extracts from HepG2 cells and 293 cells, using ADD3 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.