




AHSA1 Antibody

Product Code	CSB-PA870024
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95433
Immunogen	Synthesized peptide derived from internal of Human AHSA1.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total AHSA1 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	AHA1; activator of heat shock 90kDa protein ATPase homolog 1;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	AHSA1

Image



Western blot analysis of extracts from HT-29 cells, A549 cells and K562 cells, using AHSA1 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.